



## SYRINGE WITH CATHETER TIP

STERILE, SINGLE-USE 50 ML SYRINGE  
WITHOUT NEEDLE FOR IRRIGATION



# SYRINGE WITH CATHETER TIP

STERILE, SINGLE-USE 50 ML SYRINGE WITHOUT NEEDLE FOR IRRIGATION

- Transparent barrel for precise volume reading
- Smooth plunger motion along entire length of barrel
- Large finger grips for secure handling
- Latex-free, BPA-free, DEHP-free, PVC-free

CE<sub>0123</sub>

Class I Medical Devices  
Rule 1  
MDD 93/42/EEC  
UMDNS 13942

## PRODUCT COMPLIANCE

- CE-marked, Class I Medical Devices, Rule 1, MDD 93/42/EEC, UMDNS 13942
- Complies with the following norms, directives, and regulations:
  - ISO 7886-1
  - EN ISO 10993-1
  - EN ISO 10993-4
  - EN ISO 10993-5
  - EN ISO 10993-7
  - EN ISO 11607-1
  - EN ISO 11607-2
  - EN ISO 13485
- Labels contain 15 languages

## MANUFACTURING DETAILS

**Legal manufacturer:** Shanghai Kindly Enterprise Development Group Co. Ltd.

**Country of origin:** China

## STERILIZATION AND SHELF LIFE

### Sterilization

EtO (Ethylene oxide)



### Single-use only




### Shelf life

5 years



## MATERIALS USED



DEHP-free  
BPA-free  
Latex-free  
PVC-free

## PRODUCT RANGE OVERVIEW

Reference code	Description	Size (ml)	Type
SY3-50CTC-EC	3 part syringe 50 ml catheter tip	50-60	Centric

## PACKAGING DETAILS

Reference	Size		Quantity	Inner box			Outer box				
	L (mm)	W (mm)		L (mm)	W (mm)	H (mm)	Quantity	L (cm)	W (cm)	H (cm)	Weight (kg)
SY3-50CTC-EC	260	65	30	260	225	170	360	540	530	485	18.7

### Languages

15 languages on all levels of packaging:

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Norwegian (NO), Finnish (FI), Russian (RU), Turkish (TR), and Polish (PL)








### Transport conditions

Closed and dry

### Storage conditions

Open the packaging only immediately before use to guarantee sterility.

## Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
<b>REF</b>	✓	✓	✓
<b>LOT</b>	✓	✓	✓
	✓	✓	✓
<b>STERILE</b> <b>EO</b>	✓	✓	✓
<b>CE</b>	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
<b>EC</b> <b>REP</b>	✓	✓	✓
	—	✓	✓
	—	—	✓
	✓	✓	✓

