



Untersuchungsleuchten Examination Lamps

LED Examination Lamp

Modern practice examination lamp with energy saving LED technology. With an illuminance of approx. 50.000 lux at 30 cm working distance, the LUXAMED® examination lamp ensures optimum illumination for a professional examination.

- 15 V / 11.25 W LED Lamp with 4 high-power LEDs
- Illuminance: approx. 50.000 Lux / 30 cm working distance
- aluminium illumination head for optimal heat convection and thermal radiation
- service life: approx. 50.000 hours
- Colour temperature: approx. 4.200 K
- Gooseneck construction with 90 cm long neck, neck remains in desired position
- maximum working area: 165 cm
- robust, non-tilting 5-foot-stand, made from metal (powder-coated)
- multifunction plug (EU/UK/US)

Luxamed GmbH & Co. KG
Daniel-Weil-Straße 3
89143 Blaubeuren | Germany
T +49 (0) 7344 92905-0
F +49 (0) 7344 92905-10
info@luxamed.de
www.luxamed.de

