



Prefilled Humidifier

UWC HydroX Prefilled Humidifier is using Sterile Water. There is a diffuser designed inside the bottle to produce smaller bubbles with high humidity output to humidifying the respiratory gas for patient. Prefilled humidifier have been shown to decrease the likelihood of contamination when compared to reusable oxygen humidification units

CAUTIONS
 - Check bottle seal is intact first before use.
 - Do not use if seal is damaged



BENEFITS

- Built-in oxygen tubing connector
- Significantly increases humidity
- Gives long lasting bacteria free
- Clinical studies confirm
- Gamma sterilized
- Flow restriction alarm with humidifier adaptor

INGREDIENTS

- Sterile water
- No anti microbial agents
- No active ingredient
- No others added substances

KEMENKES RI AKL 1040 3412 156



Import & Distributor by :



Dealer :



INSTRUCTION FOR USE

